



Recombinant *Xenopus tropicalis* Protein phosphatase 1 regulatory subunit 3B (ppp1r3b)

Product Code	CSB-MP634708XBF
Abbreviation	ppp1r3b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q28E28
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MPWISSHMNS SSVLDFFPHK TMAVDIAMKF YLRSPPLRRD RVECRIARKS NEPLRPCIQT TNKTLSELS NQENKVKKRV SFADSRGLAL TMVKVYSDFD DELEIPFNIS ELIDNIVNLT TVEKERFFLD FVQPSADYLD FRNRLQAESV CLENCMLKDK ALVGTVKVKN LAFQKSVKIR MTCDSWQTYK DYDCLYVKDT YAGSDKDTFS FDVSLPEGIQ SSASIEFAVC FECEGRIFWD SNKGLNYRIV RQDHRIPSDF ESRHYDPVCL SVDQYGSPPRC SYGIFPELQT YSGFDKLGPY Y
Source	Mammalian cell
Target Names	ppp1r3b
Protein Names	Recommended name: Protein phosphatase 1 regulatory subunit 3B
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.