



Recombinant *Xenopus tropicalis* Vacuolar protein sorting-associated protein 26A (vps26a)

Product Code	CSB-BP634725XBF
Abbreviation	vps26a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q28HT6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSFLSGFFGP ICEIEVVLND ADTRKVSEIK TEEGKVEKHF LFYDGESVAG KVNIVFRQPG KRLEHQGIRI EFVGQIELFN DKSNTHEFVN LVKELALPGE LTQSRNYDFE FMQVEKPYES YIGANVRLRY FLKVTIVRRL TDLVKEYDLI VHQLATYPDV NNSIKMEVGI EDCLHIEFEY NKSXYHLKDV IVGKIYFLLV RIKIQHMELQ LIKKEITGIG PSTTTETETV AKYEIMDGAP VKGESIPIRL FLAGYDPTPT MRDVNKKFSV RYFLNLVLVD EEDRRYFKQQ EIILWRKAPE KIRKRTNFBHQ RFEPQEPQAS AEEPEI
Source	Baculovirus
Target Names	vps26a
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 26A Alternative name(s): Vesicle protein sorting 26A
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.