



Recombinant Danio rerio Eukaryotic translation initiation factor 3 subunit E-B (eif3eb)

Product Code	CSB-EP634822DIL-B
Abbreviation	eif3eb
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1LUA8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAEYDLTTRI AHFLDRHLVF PLLEFLSVKE IYNEKELLQG KLDLLSETNM VDFAMDVYKN LYPDKEIPHS LRDKRTTVVA QLKQLQSETE PIVKMFEDPE TQRQMQSTRD GRMLFEYLAD KHSFRQEYLD TLYRYAKFQY ECGNYSGAAE YLYFFRVLVP STDRNALSSL WGKLASEILM QNWEAAMEDL TRLRETIDNN TVSSPLQSLQ QRTWLIHWSL FVFFNHPKGR DNIIEFLYQ PQYLNAIQTM CPHILRYLST AVITNKDVRK RRQVLKDLVK VIQQESYTYK DPITEFVECL YVNFDFDSAQ RKLRECESVL VNDFFLVAEL EDFIENARLF IFETFCRIHQ CISISMLADK LNMTPEEAER WIVNLIRNAR LDAKIDSKLG HVVMGNNAVS PYQQVIEKTK SLSFRSQMLA MNIEKKISHS NRNETPNWAA QDTGFY
Source	E.coli
Target Names	eif3eb
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit E-B Short name= eIF3e-B Alternative name(s): Eukaryotic translation initiation factor 3 subunit 6-B
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.