



Recombinant Danio rerio TLD domain-containing protein KIAA1609 homolog (si:ch211-260p9.6, zgc:153621)

Product Code	CSB-BP634826DIL
Abbreviation	si:ch211-260p9.6, zgc:153621
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1LWV7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGNGDSVVAQ KRLARFRPDE RPAVEGTFDR LHGARSSASA GKTSKGLSLD MLQLTMGKMA SESMIKRVFQ GLHSIDPGVP LHPGDGVSRE QLLIFLADVL RGTAERAPL VLAMAEGAKA TVTTTEQIRG FMEDLVYAAV QTLAHKGHLR AWHPERMGDG AQQVKLLAEQ LTSELKPSDQ NSCDIACLED WLFRIPIPMAM FLELLIGEGL GVVLPSPPPP TLLPPCQFAP WTDLRCVLSL PLLMFLSPLL PEGHSAPWRM LFSTKMHGES FTRLLGSCKS RGPTVLLVKD TKGYIFGGFS SQSWEVKPQF QGDSRCFLFS VFPYMRVFTC TGYNDHYMYL NQQQQTMPNG LGMGGQHGYF GLWLDYDFGH GHSRARPRCT TYGSPQLSAD EDFKLDTEV WGVGKLPEEQ EEDEKKKSIL DADLEVQAMM EMTGKTLHSQ GLREPEEDED
Source	Baculovirus
Target Names	tldc1
Protein Names	Recommended name: TLD domain-containing protein KIAA1609 homolog
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.