



# Recombinant Danio rerio TLD domain-containing protein KIAA1609 homolog (si:ch211-260p9.6, zgc:153621)

<b>Product Code</b>	CSB-EP634826DIL
<b>Abbreviation</b>	si:ch211-260p9.6, zgc:153621
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1LWV7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGNGDSVVAQ KRLARFRPDE RPAVEGTFDR LHGARSSASA GKTSKGLSLD MLQLTMGKMA SESMIKRVFQ GLHSIDPGVP LHPGDGVSRE QLLIFLADVL RGTAERAPL VLAMAEGAKA TVTTTEQIRG FMEDLVYAAV QTLAHKGHRL AWHPERMGDG AQQVKLLAEQ LTSELKPSDQ NSCDIACLED WLFRIPIIMAM FLELLIGEGL GVVLPSPPPP TLLPPCQFAP WDLRCLVLSL PLLMFLSPLL PEGHSAPWRM LFSTKMHGES FTRLLGSCSKS RGPTVLLVKD TKGYIFGGFS SQSWEVKPQF QGDSRCFLFS VFPYMRVFTC TGYNDHYMYL NQQQQTMPNG LGMGGQHGYF GLWLDYDFGH GHSRARPRCT TYGSPQLSAD EDFKLDTEV WGVGKLPEEQ EEDEKKKSIL DADLEVQAMM EMTGKTLHSQ GLREPEEDED
<b>Source</b>	E.coli
<b>Target Names</b>	tldc1
<b>Protein Names</b>	Recommended name: TLD domain-containing protein KIAA1609 homolog
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.