



Recombinant Danio rerio Triosephosphate isomerase A (tpi1a)

Product Code	CSB-BP634904DIL
Abbreviation	tpi1a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1MTI4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	SSRKFFVGG NWKMNGDKES LGELIMTLNT ASLNDET DVV CGAPSIYLDY ARSKLDQRIG VAAQNCYKVP KGAF TGEISP AMIKDCGIDW VILGHSERRH VFGESDELIG QKVAHCLESD LGVIACIGEK LEEREAGTTE DVVFEQTKVI ADNVKDWTRV VLAYEPVWAI GTGKTASPEQ AQEVHEKLRG WLRANVSDAV ADSVRIIYGG SVTGGNCKEL AAQADVDGFL VGGASLKPEF VDIINARS
Source	Baculovirus
Target Names	tpi1a
Protein Names	Recommended name: Triosephosphate isomerase A Short name= TIM-A EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase A
Expression Region	2-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.