



# Recombinant Human Protein inscuteable homolog (INSC)

<b>Product Code</b>	CSB-MP634907HU
<b>Abbreviation</b>	INSC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1MX18
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRRPPGNGEA ASEGPGGWGL WGVQESRRLC CAGHDRCKQA LLQIGINMMA LPGGRHLDSV TLPGQRLHLM QVDSVQRWME DLKLMTECEC MCVLQAKPIS LEEDAQGDLI LAGGPGPGDP LQLLLKRGWV ISTELERRIGQ KLAQDRWARV HSMSVRLTCH ARSMVSEYSA VSRNSLKEMG EIEKLLMEKC SELSAVTERC LQVENEHVLC SMKACVSETL SMLGQHFQGL LELALTRETVQ ALVRKIDASD NIYTTTESTTG NLFSLTQEGA PLCRIIAKEG GVVALFKVCR QDSFRCLYPQ ALRTLASICC VEEGVHQLEK VDGVLCLADI LTDNSHSEAT RAEAAVVAQ VTSPHLPVTQ HLSSFLESME EIVTALVKLC QEASSGEVFL LASAALANIT FFDTMACEML LQLNAIRVLL EACSDKQRVD TPYTRDQIVT ILANMSVLEQ CASDIIQENG VQLIMGMLSE KPRSGTPAEV AACERVQQKA AVTLARLSRD PDVAREAVRL SCMSRLIELC RSPSERNSSD AVLVACLAAL RRLAGVCPEG LQDSDFQQLV QPRLVDSFLL CSNMEESFV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	INSC
<b>Protein Names</b>	Recommended name: Protein inscuteable homolog
<b>Expression Region</b>	1-579
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.