



# Recombinant Escherichia coli Fatty acid oxidation complex subunit alpha (fadB)

<b>Product Code</b>	CSB-BP635060EGW
<b>Abbreviation</b>	fadB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1R466
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLYKGDITYL DWLEDGIAEL VFDTPGSVNK LDTATVASLG EAIGVLEQQS DLKGLLLRSN KAAFIVGADI TEFLSLFLVP EEQLSQWLHF ANSVFNRLD LPVPTIAAVN GYALGGGCEC VLATDYRLAT PDLRIGLPET KLGIMPGFGG SVRMPRMLGA DSALEIIAAG KDVGADQALK IGLVDGVVKA EKLVEGAIAI LRQAINGDLD WKAKRQPKLE PLKLSKIEAA MSFTIAKGMV AQTAGKHYP PITAVKTIEA AARFGREEAL NLENKSFVPL AHTNEARALV GIFLNDQYVK GKAKKLTADV ETPKQAAVLG AGIMGGGIAY QSAWKGVPV MKDINDKSLT LGMTEAAKLL NKQLERKID GLKLAGVIST IHPTLDYAGF DRVDVVVEAV VENPKVKKAV LAETEQKVRP DTVLASNTST IPISLANAL ERPENFCGMH FFNPVHRMPL VEIIRGEKSS DETIAKVVAW ASKMGKTPIV VNDPCGFFVN RVLFPYFAGF SLLLRDGAGF RKIDKVMKQ FGWPMGPAYL LDVVGIDTAH HAQAVMAAGF PQRMQKDYRD AIDALFDANR FGQKNGLGFW RYKEDSKGKP KKEEDVVVDD LLAKVSQPKR DFSEEEIAR MMIPMVNEVV RCLEEGIIAT PAEADMALVY GLGFPPFHGG AFRWLDLGS AKYLDMAQQY QHLGPLYEVP EGLRNKARHN EPYYPPEVA RVPVGLKTA
<b>Source</b>	Baculovirus
<b>Target Names</b>	fadB
<b>Protein Names</b>	Recommended name: Fatty acid oxidation complex subunit alpha Including the following 2 domains: Enoyl-CoA hydratase/Delta(3)-cis-Delta(2)-trans-enoyl-CoA isomerase/3-hydroxybutyryl-CoA epimerase EC= 4.2.1.17 EC= 5.1.2.3 EC=
<b>Expression Region</b>	1-729
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.