



# Recombinant Escherichia coli Catecholate siderophore receptor Fiu (fiu)

<b>Product Code</b>	CSB-EP635133EGW
<b>Abbreviation</b>	fiu
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1REC0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LAAEGQANA DDTLVVEAST PSLYAPQQSA DPKFSRPVAD TTRTMTVISE QVIKDQGATN LTDALKNVPG VGAFFAGENG NSTTGDAIYM RGADTSNSIY IDGIRDIGSV SRDTFTEQV EVIKGPSGTD YGRSAPTGSI NMISKQPRND SGIDASASIG SAWFRRGTLT VNQVIGDTTA VRLNVMGEKT HDAGRDKVKN ERYGVAPVA FGLGTANRLY LNYLHVTQHN TPDGGIPTIG LPGYSAPSAG TAALNHSGKV DTHNFYGTDS DYDDSTTDTA TMRFEHDIND NTTIRNTTRW SRVKQDYLMT AIMGGASNIT QPTSDVNSWT WSRTANTKDV SNKILTNQTN LTSTFYTGAI GHDVSTGVEF TRETQTNYGV NPVTLPVNI YHPDSSIHG GLTRNGANAN GQTDTFAIYA FDTLQITRDF ELNGGIRLDN YHTEYDSATA CGGSGRGAIT CPAGVAKGSP VTTVDTAKSG NLVNWKAGAL YHLTENGNIY INYAVSQQPP GGNNFALAQS GSGNSANRTD FKPQKANTSE IGTKWQVLDK RLLLTAALFR TDIENEVEQN DDGTYSQYGK KRVEGYEISV AGNITPAWQM IGGYTQQKAT IKNGKDVAQD GSSSLPYTPE HAFTLWSQYQ ATDDISVGAG ARYIGSMHKG SDGAVGTPAF TEGYWVADAK LGYRVNRNLD FQLNVYNLFD TDYVASINKS GYRYHPGEPR TFLLTANMHF
<b>Source</b>	E.coli
<b>Target Names</b>	fiu
<b>Protein Names</b>	Recommended name: Catecholate siderophore receptor Fiu Alternative name(s): Ferric iron uptake protein TonB-dependent receptor Fiu
<b>Expression Region</b>	32-760
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.