



Recombinant Lactobacillus salivarius Bifunctional protein pyrR (pyrR)

Product Code	CSB-EP635240LAAM
Abbreviation	pyrR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1WTX5
Product Type	Recombinant Protein
Immunogen Species	Lactobacillus salivarius (strain UCC118)
Purity	≥85% (SDS-PAGE)
Sequence	MAKQIYDSMA MKRALTRMTY EIIKKNKGID DLVLVGKTR GVYLAQRIAD RLKQLENAEV PVGQLDITLY RDDRHDA SLA QDPVVNEADL GFDIKDKHVI LVDDVLFTGR TIRAALDALM DIGRPKKINL AVLVDRGHRE LPIRADFVGK NIPTAMDEQV AVYVDEV D GK DGIELKKL
Source	E.coli
Target Names	pyrR
Protein Names	Recommended name: Bifunctional protein pyrR Including the following 2 domains: Pyrimidine operon regulatory protein Uracil phosphoribosyltransferase Short name= UPRTase EC= 2.4.2.9
Expression Region	1-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.