



Recombinant Fatty acid-binding protein homolog 2 (lbp-2)

Product Code	CSB-BP635304CXY
Abbreviation	lbp-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20224
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	E QLPEKFGYGTFLDLHSENFDE YLTAKGYGWF TRKLVTFATF KKVFANKANK NLFDYSNLTS KKDVFYKNVQ IGSKFEGEGL DNTKHEVTFT LKDGHLFEHH KPLEEGESKE ETYEYYFDGD FLIQKMSFNN IEGRRFYKRL P
Source	Baculovirus
Target Names	lbp-2
Protein Names	Recommended name: Fatty acid-binding protein homolog 2 Alternative name(s): Lipid-binding protein 2
Expression Region	20-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.