



Recombinant Bovine Angiopoietin-4 (ANGPT4)

Product Code	CSB-EP635541BO
Abbreviation	ANGPT4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24K15
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	QRRGQEEA GRRRAHRVQH GQCSYTFVLP EPEPCPEPE AFGGSNSLQR DSPAATLNLG DWPSQRMRL EKMLENNTQW LQKLERIYIQV NLRLELAQAQ QHMQNQQTAT MLELGTSLLT QTTAQTRKLT DVEAQVLNQT SRMEIQLLET SLSTNKLEKQ LLLQGHELHR LQGHNSALET RVQALETQQQ AELASLSGEK ERLRLLGRQ SGALAGLERT LRAASSNSSL LQRQQHQLLE SVQRLVRVMA QGPASMRAAD QLFQDCAEIQ RFGANASGIY TIHVANVTEP RKFVCDMEAS GGGWTLIQR ENGSVNFQRN WKDYKQGFGN PAGEHWLGNE VVHQLTSRAT YSLRVELQDW EGNEAYAQYE HFQLGSEAQL YRLSLSGYSG SAGRQSSLVL QGTNFFSTRDA DNDNCLCKCA QMLSGGWWFD ACGLSNLNGI YYPARHHVRK LNGIRWHYFQ GPSYSLRTRR MMVRPSGI
Source	E.coli
Target Names	ANGPT4
Protein Names	Recommended name: Angiopoietin-4 Short name= ANG-4
Expression Region	23-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The mechanism by which they contribute to angiogenesis is thought to involve regulation of endothelial cell interactions with supporting perivascular cells. This protein functions as an agonist and is an angiopoietin.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.