



Recombinant Enolase (enol-1)

Product Code	CSB-MP635655CXY
Abbreviation	enol-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27527
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	PITKIHARQ IYDSRGNPTV EVDLFTEKGV FRAAVPSGAS TGVHEALELR DGDKAVHLGK GVLKAVSNIN EKIAPALIAK GFDVTAQKDI DDFMMALDGS ENKGNL GANA ILGVSLAVAK AGAVHKGLPL KYIAELAGT GKVVLPVPAF NVINGGSHAG NKLAMQEFMI LPVGASSFAE AMRMGSEVYH HLKAEIKKRY GLDATAVGDE GGFAPNIQDN KEGLDLLNTA IDKAGYTGKI SIGMDVAASE FFKDGKYDLD FKNPASDSSK WLSGEQLTEL YQSFIKEYPV VSIEDAFDQD DWDNWKFGHG ATSIQLVGDD LTVTNPKRIQ TAIDKKSCNC LLLKVNQIGS VTESIEAAKL SRANGWGMV SHRSGETEDT FIADLVVGLA TGQIKTGAPC RSERLAKYNQ LLRIEEELGA DAVYAGHNFR NPQV
Source	Mammalian cell
Target Names	enol-1
Protein Names	Recommended name: Enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2-phosphoglycerate dehydratase
Expression Region	2-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.