



Recombinant Chicken Annexin A2 (ANXA2)

Product Code	CSB-BP001840CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17785
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	STVHEILSK LSLEGDHSLP PSAYATVKAY SNFDADRDA ALEAAIKTKG VDEVTIINIL TNRSNEQRQD IAFAYQRRTK KELSAALKSA LSGHLEAVIL GLLKTPSQYD ASELKAAMKG LGTDEDTLIE IICSR TNQEL NEINRVYREM YKTELEKDII SDTSGDFRKL MVALAKGKRC EDTSV IDYEL IDQDARELYD AGVKRKGTDV PKWINIMTER SVP HLQKVFE RYKSYSPYDM LESIKKEVKG DLENAFLNLV QCIQNKQLYF ADRLYDSMKG KGTRDKVLIR IMVSRCEVDM LKIKSEFKRK YGKSLYYFIQ QDTKGDYQRA LLNLGGED
Source	Baculovirus
Target Names	ANXA2
Protein Names	Recommended name: Annexin A2 Alternative name(s): Annexin II Annexin-2 Calpactin I heavy chain Calpactin-1 heavy chain Chromobindin-8 Lipocortin II Placental anticoagulant protein IV Short name= PAP-IV Protein I
Expression Region	2-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.