



# Recombinant Chicken Annexin A2 (ANXA2)

<b>Product Code</b>	CSB-YP001840CH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17785
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	STVHEILSK LSLEGDHSLP PSAYATVKAY SNFDADRDA ALEAAIKTKG VDEVTIINIL TNRSNEQRQD IAFAYQRRRK KLSAALKSA LSGHLEAVIL GLLKTPSQYD ASELKAAMKG LGTDEDTLIE IICSR TNQEL NEINRVYREM YKTELEKDII SDTSGDFRKL MVALAKGKRC EDTSVI DYEL IDQDARELYD AGVKRKGTDV PKWINIMTER SVPHLQKVFE RYKSYSPYDM LESIKKEVKG DLENAFLNLV QCIQNKQLYF ADRLYDSMKG KGTRDKVLIR IMVSRCEVDM LKIKSEFKRK YGKSLYYFIQ QDTKGDYQRA LLNLGGED
<b>Source</b>	Yeast
<b>Target Names</b>	ANXA2
<b>Protein Names</b>	Recommended name: Annexin A2 Alternative name(s): Annexin II Annexin-2 Calpactin I heavy chain Calpactin-1 heavy chain Chromobindin-8 Lipocortin II Placental anticoagulant protein IV Short name= PAP-IV Protein I
<b>Expression Region</b>	2-339
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.