



Recombinant *Rhizobium leguminosarum* bv. *viciae* Lipoyl synthase (lipA)

Product Code	CSB-EP635883RKS
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1MH26
Product Type	Recombinant Protein
Immunogen Species	<i>Rhizobium leguminosarum</i> bv. <i>viciae</i> (strain 3841)
Purity	>85% (SDS-PAGE)
Sequence	MVTILDITNP DAKRVRHPEK AHRPDTEVMR KPDWIRVKAP TSKGYAETRA IVKEHKLTV CEEAGCPNIG ECWDKKHATF MIMGEICTRA CAFCNVATGK PNALDMAEPE NVAKAVKEMG LSHVVITSVD RDDLEDGGAE HFEKVIWAIR AASPMTTIEI LTPDFLKKPG ALERVVAAKP DVFNHNMETV PGNYLTVRPG ARYFHSVRLQ QRVKELDPTM FTKSGIMVGL GEERNEVLQL MDDLRTADVD FLTIGQYLQP TRKHHKVESF VTPDEFKSYE TVAYSKGFLM VASSPLTRSS HHAGDDFARL RAAREKLLM AAE
Source	<i>E.coli</i>
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.