



Recombinant Danio rerio Coiled-coil domain-containing protein 104 (ccdc104)

Product Code	CSB-MP636236DIL
Abbreviation	ccdc104
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1RM35
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAEDSEWVPE SIAGFLGSPD WLIPLAEFME SRCSVFDDDED ENKLTYTEIH QQYKQLVERL LENYMQDVGI SEQQFLEACS SFSSSASLQT VFQPVLATDD FQMRSLMVQ KNMELQLQAL HVIKQRNGGL PECLTDGVDV MSELEQQEMK ILQEVLRRSK EYDLQMGQC GLGSEDLAST SSSVSETPQN PEQRLSNGVS DPLTLTQPDS EMEESSTATQ RKMEESEPTA TACKSKALP AVRAPLRGSV VAKQSSAERE RELERARDTH TDCRTLTDDEV ALQQRSEYLK QQRDKLQALK RQTGKGLTDT AAAAPEPKP ATQEISVEEK KKLQKRKHVA EKLKEEVIKK
Source	Mammalian cell
Target Names	cfap36
Protein Names	Recommended name: Coiled-coil domain-containing protein 104
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.