



# Recombinant Bovine Ubiquitin-conjugating enzyme E2 S (UBE2S)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP636238BO  |
| <b>Abbreviation</b>      | UBE2S   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q1RML1  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Bos taurus (Bovine)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MNSNVENLPP HIIRLVYKEV TTLTADPPDG IKVFPNEEDL TDLQVTIEGP<br>EGTPYAGGLF RMKLLL GKDF PASPPKGYFL TKIFHPNVGA NGEICVNVLK<br>RDWTAELGIR HVLLTIKCLL IHPNPESALN EEAGRLLLEN YEEYAARARL<br>LTEIHGGAGG PSGGRPEPGR ATASGAAAST ADPTAPGGPA GAEGPMAKKH<br>AGERDKKLAA KKKTDKKRAL RRL  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | UBE2S   |
| <b>Protein Names</b>     | Recommended name: Ubiquitin-conjugating enzyme E2 S EC= 6.3.2.19<br>Alternative name(s): Ubiquitin carrier protein S Ubiquitin-protein ligase S   |
| <b>Expression Region</b> | 1-223   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.