



Recombinant 26S proteasome non-ATPase regulatory subunit 6 (rpn-7)

Product Code	CSB-EP636336CXY
Abbreviation	rpn-7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20585
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTEAAKKSTK KPVDDGNFDK EIISRWPDL ELSQTRFMLNH PEVDSSVKEA KLEKLQETIK EFDMAPFYEL VCADFKIVVD ATQLAAMKAA NQKKIDEITA EVEDAEKNLG ESEVRQGLLR KFEYYCQIGD KDNALKAYTA TYEKTVGMGY RIDVVFAMIR VGLFFLDHHL INKFITKAKE LMEQGGDWER KNRLRSYEAL YRMSVRDFAG AADLFLEAVP TFGSYELMTY ENLILYTVIT TTFALDRPDL RTKVIRCNEV QEQLTGGGLN GTLIPVREYL ESYYDCHYDR FFIQLAALES ERFKFDRLYS PHFNYYSRGM RHRAYEQFLT PYKTVRIDMM AKDFGVSRF IDRELHRLIA TGQLQCRIDA VNGVIEVNHR DSKNHLYKAV IKDGDILLNR IQKLARVINA
Source	E.coli
Target Names	rpn-7
Protein Names	Recommended name: 26S proteasome non-ATPase regulatory subunit 6 Alternative name(s): 26S proteasome regulatory subunit rpn-7
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.