



Recombinant Putative pyridoxamine 5'-phosphate oxidase (F57B9.1)

Product Code	CSB-EP636342CXY-B
Abbreviation	F57B9.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20939
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	METPSIDIQN IRAKYLNSHD PYLLESKLPT TSPFELFDIW FRNVASQS DLTFEEINAVSL STVGKDLRPS SRMVLLKAYT PTGFSFYTNY TSRKGNQLEE NPNAAMLFYW PKVNRQIRVE GVVEKLPDEM AVAYWNSRPV ASRIGSKSSD QSKVVPDREF LESKKVALTE LSVREGAQAI TKPESWGGYH LIPRYFEFWQ GQSDRLHDRI VFERDVDVWL LKRLSP
Source	E.coli
Target Names	F57B9.1
Protein Names	Recommended name: Putative pyridoxamine 5'-phosphate oxidase EC= 1.4.3.5 Alternative name(s): PNP/PMP oxidase Short name= PNPOx
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.