



Recombinant Uncharacterized protein R102.4 (R102.4)

Product Code	CSB-MP636399CXY
Abbreviation	R102.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21890
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MISRHFTVLT ALNLQTKIKI CIRNFSIRGA MTKVNGNSNG NGISASTDLL DTPSAYTQKS NKTHTSIDLR SDTVTVPSVE MRRAMAEIV GDDVYGEDTT TNRLEQRCAE LFGKEAGLFV TSGTMGNLLA IMAHCQRGEE IIVGRYNHIH RWEQGNYAQF AGISATTLEV KPDGTMDLND IEQAIRVKDC HMPASKLICI ENTHNYTGGK ALPIEWMRSV KQLAERRDLK VHMDGARIYN AAVASNCSVS KIASFADTVQ MCFSKGLGAP VGSIVVGPDK FIDRARHSRK ALGGGWQRQSG ILAAAHAIAL DHADATIRAD HERAKTLARM INDATPEEFR TKVFAAEKDI TNMVLVHCQN GVTVQQLTDF FQKHDILAMT FDARRIRMLV NWNVSDENLE TIVEVYKKFL KQL
Source	Mammalian cell
Target Names	R102.4
Protein Names	Recommended name: Uncharacterized protein R102.4
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.