



Recombinant Uncharacterized protein T19C3.1 (T19C3.1), partial

Product Code	CSB-MP636519CXY
Abbreviation	T19C3.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22566
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MYIANYFIQI HRFLCFSQIL NEMTVSYQLD VSSGNPLLFL RLLGRWRGSI WKSVVGDLFV WLLFYAIYF AYRYAFSKQL QTVFEEISIH TDDRMKYLPL TFMLGFFVTT VFERWRSALN VMPFIESVAL SVAVLLPGKG REDRLTRRAI IRYVVLHQIL VFRDISMRVR RRFPTLKYVV DAGFMRQEEL DVLESVNQES SQTYWVPINW ANSLALVAHQ QKLIDQPTAF NNVIFAIKEF RVAMETLIK DAIPIPIAYP QVVFLAVRVY FAICLVSRQF LISDMKSKTQ MDWPVPIMTV LEFIFVIGWM KVAEVLLNPL GEDDDDFEVN SIIDNNISRG MAIVDTTHGY HPDLVDDVFS DPNYLPAYSE NSQIPRNLTG SAAKVELAAP TDEVKIVRVN PEDAPPTSER GSMFTRRNAY SLRNRKISIA SNNLESPTSQE RKFNL SMPAG MLNKSTQPDR PTMETVSEEH EPSHFYRGDR VHSSDSVIPS FRRSFFQIKA F
Source	Mammalian cell
Target Names	best-17
Protein Names	Recommended name: Uncharacterized protein T19C3.1
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.