



Recombinant BRO1 domain-containing protein BROX homolog (B0507.2)

Product Code	CSB-EP636522CXY-B
Abbreviation	B0507.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22885
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSHWFHRNPI KPTEFVKFDL KGVLTDTCS KICGELRLRR DKLVSQFKNA SNDLEEVTK E FNEYLR LFAG FLIEIQSSMV ELENKDAGNK NSKLIPLIRF KWGNSMLPQA ATEVS DTWFE ALSMIQCMAM WLTKHAASMA GKDEVRESDA KECLQCLRQA GG MFQYVKDE SSRLSGANEV EGSDFDPKVM ETYILTATAE AQEVIVARAI EMKHDDGLIS SLAAVTASIF SKADQCLNNL PDES FARWRR YLQLKHHFY L AYAF A FLGQK QLS EDKCGEA VRACKQGIAE YGVAKEMAAM YATATGPGTR IKPEQHLFFR RIEPLLNRHL EKAERENGFI YHQKVPDEIP QLDVEATYGL AKLDSFTYPP PAESWNTAVY SAFDLSKANM PDFSKIKKSK SKLDPVHEEK IYQTEKDPSN SSGCVIA
Source	E.coli
Target Names	B0507.2
Protein Names	Recommended name: BRO1 domain-containing protein BROX homolog
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.