



Recombinant Platelet-activating factor acetylhydrolase homolog 2 (paf-2)

Product Code	CSB-EP636523CXY
Abbreviation	paf-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22943
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGSYISSPQV LTRQVSGQFQ VGCKDLMIDG TVLGDRGLFM RLYFPTDSQA ADISSYPLWL PKPQYAHGLG EYLGQSSQKM NVITSTVVGE KREDCIENAQ MSTKCDKWPI VVFSHGLGGS RTFYSTYCTS LASHGYVVA VEHKDHSACW TYQLTEKNGE LVEQPIKIKL IEKNEKNEFK IRNQVVGKRV TECVKALNVL EQLNLGTVPE KVLIGNDYNW AQFKNKLVMS SASVIGHSFG GATSLASSAY TTDFQKAIVF DGWMYPLDST QQEQAKQPTL FLNVGDWQWN ENLDVMKKII SHNDGNLALT LNGAVHQCFS DFPFIFPSWL AKKFGVQGRT EPSLCMQAAI ELSLAFLENG KDGAQKCLKDE KFSSFISNEI YGREKYKL
Source	E.coli
Target Names	paf-2
Protein Names	Recommended name: Platelet-activating factor acetylhydrolase homolog 2 EC=3.1.1.47
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.