



Recombinant Drosophila melanogaster Splicing factor U2AF 50 kDa subunit (U2af50)

Product Code	CSB-BP636558DLU
Abbreviation	U2af50
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24562
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGYDDRERDR ERRRHRSRSR DRHRERSRDR RHRNSRRKP SLYWDVPPPG FEHITPMQYK AMQASGQIPA SVVPDTPQTA VPVVGSTITR QARRLYVGNI PFGVTEEMM EFFNQQMHLV GLAQAAGSPV LACQINLDKN FAFLEFRSID ETTQAMAFDG INLKGQSLKI RRP HDYQPMP GITDTPAIKP AVVSSGVIST VVPDSPHKIF IGGLPNYLND DQVKELLSF GKLRAFNLVK DAATGLSKGY AFCEYVDLSI TDQSIAGLNG MQLGDKKLIV QRASVGAKNA QNAANTTQSV MLQVPGLSNV VTSGPPTTEVL CLLNMVTPDE LRDEEYEDI LEDIKEECTK YGVVRSVEIP RPIEGVEVPG CGKVFVEFNS VLDCQKAQQA LTGRKFSDRV VVTSYFDPDK YHREF
Source	Baculovirus
Target Names	U2af50
Protein Names	Recommended name: Splicing factor U2AF 50 kDa subunit Alternative name(s): U2 auxiliary factor 50 kDa subunit U2 snRNP auxiliary factor large subunit
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.