



Recombinant *Desulfitobacterium hafniense* Enolase 1 (eno1)

Product Code	CSB-MP636598DAAE
Abbreviation	eno1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24YW4
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfitobacterium hafniense</i> (strain Y51)
Purity	>85% (SDS-PAGE)
Sequence	MEGFKIKSVK ARQVFDSRAN PTVEVDITLK DGTVGRGIVP SGASTGLFEA VELRDGEETF DGKGVSKAIH HVNEQLAPLL LDMDVRDQKA IDSRLIEIDG TPNKSRLGAN AILGCSMAAC WAGANYKVP LYRYLGGSAA NKLPVPMVQI IGGGAHADNV IDIQDFLVIP TSAPTFSEGY EMVVNVYNA KMIKFGAGKP VSIADGGGLW PTGFQSNEEG LKLLCESIEL AGYTPGQDLG IALDIASSEF YDPERGTYRL ELEKKTLTKE EMVGMLSDWV DNYPIISIED GMSELDWEGN LLLTQKLGKK IQLIGDDLFT TNIERIRKGV EMKVDNAVLI KMNQIGTITE TITIEFTQN HGYPVVSAR SGETEDCTIV HLAIATNAGQ LKVGSAARSE RTAKWNEVLR IEEALGSSAR YPSKGVFAAA GIQFNQY
Source	Mammalian cell
Target Names	eno1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.