



Recombinant V-type proton ATPase subunit B (VAPB)

| | |
|--------------------------|--|
| Product Code | CSB-EP636646PLO-B |
| Abbreviation | VAPB |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q25691 |
| Product Type | Recombinant Protein |
| Immunogen Species | Plasmodium falciparum |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSKEVVNTKA EASRVNALAA VRNYKVCPR L EYKTISGVQG PLVIIEDVKF PKYSEIVTIIH LSDNTTRQGQ ILEVCGKKAV IQVFEGTSGI DNKNSYVEVS GDILKMPMSD EMLGRVFNGS GKPIDKGPNI LADDYLDING NPINPQCRVY PKEMIQTGIS TIDVMNSIVR GQKIPLFSAA GLPHNEIGAQ ICRQASLVQG KDVLHDHSDDN FAVVFGAMGV NMETARYFRQ DFEENGKMER VCLFLNLAND PTIERILTPR IALTTAEYLA FEKEMHVFVI LTDMSSYADA LREVSSAREE VPGRRGYPGY MYSDLSTIYE RAGRVEGRNG SITQFPILTM PNDDITHPIP DLTGYTEGQ IFVDRNLYNR QIYPPINVLP SLSRLMKSGI GHNMTRIDHP YVSDQLYSNY AIAQDVKAMK AVIGEEALSN DDILYLEFLD KFEKRFITQN TYECRDIYQS LDIAWELLRI FPEDMLKKIK TDILSKYYPR HHAN |
| Source | E.coli |
| Target Names | VAPB |
| Protein Names | Recommended name: V-type proton ATPase subunit B Short name= V-ATPase subunit B Alternative name(s): Vacuolar proton pump subunit B |
| Expression Region | 1-494 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.