



Recombinant V-type proton ATPase subunit B (VAPB)

Product Code	CSB-EP636646PLO
Abbreviation	VAPB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q25691
Product Type	Recombinant Protein
Immunogen Species	Plasmodium falciparum
Purity	>85% (SDS-PAGE)
Sequence	MSKEVVNTKA EASRVNALAA VRNYKVCPR L EYKTISGVQG PLVIIEDVKF PKYSEIVTIIH LSDNTTRQGQ ILEVCGKKAV IQVFEGTSGI DNKNSYVEVS GDILKMPMSD EMLGRVFNGS GKPIDKGPNI LADDYLDING NPINPQCRVY PKEMIQTGIS TIDVMNSIVR GQKIPLFSAA GLPHNEIGAQ ICRQASLVQG KDVLDHSDDN FAVVFGAMGV NMETARYFRQ DFEENGKMER VCLFLNLAND PTIERILTPR IALTTAEYLA FEKEMHVFVI LTDMSSYADA LREVSSAREE VPGRRGYPGY MYSDLSTIYE RAGRVEGRNG SITQFPILTM PNDDITHPIP DLTGYTEGQ IFVDRNLYNR QIYPPINVLP SLSRLMKSGI GHNMTRIDHP YVSDQLYSNY AIAQDVKAMK AVIGEEALSN DDILYLEFLD KFEKRFITQN TYECRDIYQS LDIAWELLRI FPEDMLKKIK TDILSKYYPR HHAN
Source	E.coli
Target Names	VAPB
Protein Names	Recommended name: V-type proton ATPase subunit B Short name= V-ATPase subunit B Alternative name(s): Vacuolar proton pump subunit B
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.