



# Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At4g33900 (At4g33900)

<b>Product Code</b>	CSB-MP636977DOA
<b>Abbreviation</b>	At4g33900
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1PE27
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEYGEEPSIK RFLMLPDDL V FNCLARVSRL HYPTLSLVSK KFRFLLASKE LYQTRILLGG TESCLYVCVR LHTDSEQLHW FIIYQGNSS KKVLPPISSP NFTSAALPGF VVVGHEIYAI GGGSENKNAS INATGSKTYN ALSSVMVMDS RSHTWREAPS MRVARVFPSA CTLDGRIYVT GGCENLN SMN WMEIFDTKTQ TWEFLQIPSE EVCKGSEYLS ISYQRTVYVG SREKDVTYKM HKGKWRGADI CLNHGWSLDP SSCCVIENVF YRCSLGDVRW YDLKKREWAA LKGLEGLPTF TNYRNFKSA DHCGLAISW EEYVLVDEET KIWCAEIAIQ KRQNGEIWGT LEWFDNVFIS SGPNRHVDLL VNALTATVW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At4g33900
<b>Protein Names</b>	Recommended name: F-box/kelch-repeat protein At4g33900
<b>Expression Region</b>	1-379
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.