



# Recombinant Escherichia coli 3-octaprenyl-4-hydroxybenzoate carboxy-lyase (ubiD)

<b>Product Code</b>	CSB-BP637128EGW
<b>Abbreviation</b>	ubiD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1R469
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDAMKYNDLR DFLTLEQQG ELKRITLPVD PHLEITEIAD RTLRAGGPAL LFENPKGYSM PVLCLNLFGTP KRVAMGMGQE DVSALREVGK LLAFLKEPEP PKGFRDLFDK LPQFKQVLNM PTKRLRGAPC QQKIVSGDDV DLNRIPIMTC WPEDAAPLIT WGLTVTRGPH KERQNLGIYR QQLIGKNKLI MRWLSHRGGA LDYQEWCAAH PGERFPVSA LGADPATILG AVTPVPDTLS EYAFAGLLRG TKTEVVKCIS NDLEVPASAE IVLEGYIEQG ETAPEGPYGD HTGYINEVDS FPVFTVTHIT QREDAIYHST YTGRPPDEPA VLGVALNEVF VPILQKQFPE IVDFYLPPEG CSYRLAVVTI KKQYAGHAKR VMMGVWSFLR QFMYTKFVIV CDDVDNARDW NDVIWAITTR MDPARDTVLV ENTPIDYLDL ASPVSGLGSK MGLDATNKWP GETQREWGRP IKKDPDVAH IDAIWDELAI FNNGKSA
<b>Source</b>	Baculovirus
<b>Target Names</b>	ubiD
<b>Protein Names</b>	Recommended name: 3-octaprenyl-4-hydroxybenzoate carboxy-lyase EC=4.1.1.- Alternative name(s): Polyprenyl p-hydroxybenzoate decarboxylase
<b>Expression Region</b>	1-497
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.