



Recombinant Escherichia coli UPF0042 nucleotide-binding protein YhbJ (yhbJ)

Product Code	CSB-MP637149EGW
Abbreviation	yhbJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1R6C9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MVLMIVSGRS GSGKSVALRA LEDMGFYCVD NLPVVLLPDL ARTLADREIS AAVSIDVRNM PESPEIFEQA MSNLPDAFSP QLLFLDADRNL TLRYSRSDTR RLHPLSSKNL SLESAIDKES DLLEPLRSRA DLIVDTSEMS VHELAEMLRT RLLGKREREL TMVFESFGFK HGIPIDADYV FDVRFNPNH WDPKLRPMTG LDKPVA AFLD RHTEVHNFYIY QTRSYLELWL PMLETNNSY LTVAIGCTGG KHRSVYIAEQ LADYFRSRGK NVQSRHRTLE KRKP
Source	Mammalian cell
Target Names	rapZ
Protein Names	Recommended name: UPF0042 nucleotide-binding protein YhbJ
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.