



Recombinant Escherichia coli Crotonobetainyl-CoA:carnitine CoA-transferase (caiB)

Product Code	CSB-MP637216EGW
Abbreviation	caiB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1RGG0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDHLPMPKFG PLAGLRVVFS GIEIAGPFAG QMFAEWGAEV IWIENVAWAD TIRVQPNYPQ LSRRNLHALS LNIFKDEGRE AFLKLMETTD IFIEASKGPA FARRGITDEV LWQHNPCLVI AHLSGFGQYG TEEYTNLPAY NTIAQAFSGY LIQNGDVDQP MPAFPYTADY FSGLTATTAA LAALHKVRET GKGESIDIAM YEVMLRMGQY FMMDYFNGGE MCPRMTEGKD PYYAGCGLYK CADGYVMEL VGITQIAECF KDIGLAHLLG TPEIPEGTQL IHRIECPYGP LVEEKLDLAWL AAHTIAEVKE RFAELNIACA KVLTVPELES NPQYVARES TQWQTMGRT CKGPNIMPKF KNNPGQIWRG MP SHGMDTAA ILKNIGYSEN DIQELVSKGL AKVED
Source	Mammalian cell
Target Names	caiB
Protein Names	Recommended name: Crotonobetainyl-CoA:carnitine CoA-transferase EC=2.8.3.-
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.