



Recombinant Human Protein FAM110C (FAM110C)

Product Code	CSB-MP637287HU
Abbreviation	FAM110C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1W6H9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRALAALSAP PNERLLPRDP AATRDPDAAR PARRSAVERL AADRAKYVVRG RPGTGRGVAS EGSGPGAIC PGNDPGPPAR APAPVARRAI ARKPLRPDSL IYRQKCEFV RGSGADGPRA SLVKKLFQGP GKDKAPVPRT GDEGKAGNPE TVPTTPGPAA DPAIPETPAP AARSAAPSSV PAAPPGPEPR VVRRRRLQRS QSDLSSRYSA ALAESDTFFQ YCGLDPEVVE ALGRENFTAG SDCVTLKVRS VSVATSGSGF SRHSGGDDEG LQEEELIEQV PSTTSVIERN ARIKWLYTC KKAKETPSQE QSRTRGSKPS R
Source	Mammalian cell
Target Names	FAM110C
Protein Names	Recommended name: Protein FAM110C
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.