



Recombinant MAP kinase kinase mkk-4 (mkk-4)

Product Code	CSB-BP637374CXY
Abbreviation	mkk-4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20347
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MVQEDDENLR NSMSLRPTSL STRPTSLSVN GNEKTLPEES VLRSLSSTGTL KYPDDEHLYT FSSANLQDLG AIGNGNFGTV YKMRHKETGK LIAVKRIRCN NIGHREQIRL LREHDTIVKS EKGPNIVKFY GAIFSEGDCW ICMELMDISM DLLYKRVYMV KNSRLNENVV GHITVCTVDA LDYLKKEKELKI IHRDVKPSNI LVDGTGAVKL CDFGICGQLE ESFAKTHDAG CQPYLAPERI TSSDKYDVRS DVWSLGITLY EIATGKFPYQ EWNSLFDQIA TVVSGDPPIL HPDSDDFHYS LPLVKFINTC LTKDRRHRPK YDTLKSDFDY RIYAVAGPEI EEAKRILGVE AIDTRNHPVD HRG
Source	Baculovirus
Target Names	mkk-4
Protein Names	Recommended name: MAP kinase kinase mkk-4 EC= 2.7.12.2
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.