



Recombinant Probable 39S ribosomal protein L49, mitochondrial (R11D1.9)

Product Code	CSB-MP637440CXY
Abbreviation	R11D1.9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21939
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLGSSLRIFS PRAVSRLSFS TTSVEKAEKL WENPWKHALP EKSRSLATVE EAPIDWSYVE RLMPIEVVPN VPEHEKYPTP SGWTPPTEAA KTHQYYIRRR HDHLLPLYLE RKRDLLNEKT LDFDYVELVT IRNVDGDIFA CENDLRSYLE EHLGHSIASH VDELKGRIKI KGAPRVLIEQ FFYSKGF
Source	Mammalian cell
Target Names	mrpl-49
Protein Names	Recommended name: Probable 39S ribosomal protein L49, mitochondrial Short name= L49mt Short name= MRP-L49
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.