



Recombinant *Drosophila melanogaster* Dynein light chain 1, cytoplasmic (ctp)

| | |
|--------------------------|---|
| Product Code | CSB-EP637582DLU-B |
| Abbreviation | ctp |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q24117 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila melanogaster</i> (Fruit fly) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MSDRKAVIKN ADMSEEMQQD AVDCATQALE KYNIEKDIAA YIKKEFDKKY NPTWHCIVGR NFGSYVTHET RHFYFYLGQ VAILLFKSG |
| Source | E.coli |
| Target Names | ctp |
| Protein Names | Recommended name: Dynein light chain 1, cytoplasmic Alternative name(s): 8 kDa dynein light chain Cut up protein |
| Expression Region | 1-89 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |