



Recombinant *Drosophila melanogaster* Dynein light chain 1, cytoplasmic (ctp)

Product Code	CSB-EP637582DLU
Abbreviation	ctp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24117
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSDRKAVIKN ADMSEEMQQD AVDCATQALE KYNIEKDIAA YIKKEFDKKY NPTWHCIVGR NFGSYVTHET RHFYFYLGQ VAILLFKSG
Source	<i>E.coli</i>
Target Names	ctp
Protein Names	Recommended name: Dynein light chain 1, cytoplasmic Alternative name(s): 8 kDa dynein light chain Cut up protein
Expression Region	1-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.