



Recombinant *Drosophila virilis* Hsp90 co-chaperone Cdc37 (Cdc37)

Product Code	CSB-YP637600DMP
Abbreviation	Cdc37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24740
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVDYSKWKNI EISDDEDETH PNIDTPSLFR WRHQARVERM AESEKEKEEL KKKRQSYQAR LINVKELISK KEGDEVALKK ELEKIENEGK ELDRQENELL KREKTPWNV DTISKPGFEK TVINKKPPRK AAEDLTEER EQRMKQFVKE NEKLCKQYGM LRKYDDSKRF LQEHLDLVCE ETANYLVIWS INLEMEEKHD LMAHVAHQCI CMQYILELAK QLDVDPRACV SSFFSKIQQC VPEYRQQFES EIEGFKERIQ KRAQEKLQEA MAQLEEEEEKQ ERMGPGGLDP ADVFESLPEE LKACFESRDI ELLQKTIATM PVDQAKYHMK RCVDSGLWVP NAADVQLGND ETEAKDDETT TATNETPTEK ENSKSQNEKS ADAAKEEPLY AGVSTDDVD
Source	Yeast
Target Names	Cdc37
Protein Names	Recommended name: Hsp90 co-chaperone Cdc37 Alternative name(s): Hsp90 chaperone protein kinase-targeting subunit
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.