



Recombinant *Drosophila melanogaster* Protein deadpan (dpn)

Product Code	CSB-BP637696DLU
Abbreviation	dpn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q26263
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDYKNDINSD DDFDCSNGYS DSYGSNGRMS NPNGLSKAEL RKTNKPIMEK RRRARINHCL NELKSLILEA MKKDPARHTK LEKADILEMT VKHLQSVQRQ QLNMAIQSDP SVVQKFKTGF VECAEEVNRY VSQMDGIDTG VRQLSAHLN QCANSLEQIG SMSNFSNGYR GGLFPATAVT AAPTPLFPSL PQDLNNSRT ESSAPAIQMG GLQLIPSRP SGEFALIMPN TGSAAPPPGP FAWPGSAAGV AAGTASAALA SIANPHTLND YTQSFMSAF SKPVNTSVPA NLPENLIHTL PGQTQLPVKN STSPPLSPIS SISSHCEESR AASPTVDVMS KHSFAGVFST PPPTSAETSF NTSGSLNLSA GSHDSSGCSR PLAHLQQQV SSTSGLAKRD REAEAESSDC SLDEPSSKKF LAGAIEKSSS AWRPW
Source	Baculovirus
Target Names	dpn
Protein Names	Recommended name: Protein deadpan
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.