



# Recombinant Drosophila melanogaster Protein deadpan (dpn)

<b>Product Code</b>	CSB-EP637696DLU-B
<b>Abbreviation</b>	dpn
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q26263
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDYKNDINSD DDFDCSNGYS DSYGSNGRMS NPNGLSKAEL RKTNKPIMEK RRRARINHCL NELKSLILEA MKKDPARHTK LEKADILEMT VKHLQSVQRQ QLNMAIQSDP SVVQKFKTGF VECAEEVNRY VSQMDGIDTG VRQLSAHLN QCANSLEQIG SMSNFSNGYR GGLFPATAVT AAPTPLFPSL PQDLNNSRT ESSAPAIQMG GLQLIPSRP SGEFALIMPN TGSAAPPPGP FAWPGSAAGV AAGTASAALA SIANPHTLND YTQSFMSAF SKPVNTSVPA NLPENLIHTL PGQTQLPVKN STSPPLSPIS SISSHCEESR AASPTVDVMS KHSFAGVFST PPPTSAETSF NTSGSLNLSA GSHDSSGCSR PLAHLQQQV SSTSGLAKRD REAEAESSDC SLDEPSSKKF LAGAIEKSSS AWRPW
<b>Source</b>	E.coli
<b>Target Names</b>	dpn
<b>Protein Names</b>	Recommended name: Protein deadpan
<b>Expression Region</b>	1-435
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.