



Recombinant *Drosophila melanogaster* Proteasome subunit alpha type-7-1B (Pros28.1B)

Product Code	CSB-EP637725DLU-B
Abbreviation	Pros28.1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27575
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAQRYDRAVT IYSPDGHLLQ VEYAQEAVRR GSTVMGLRTN NAIVIGVEKR SVGDLQEERM VRKICMLDDH VVMTFSGTLTA DARILVSRAQ MEAQSHRLNF EKPTTVEYIT RYIAQLKQNY TQSNRRPFG LSCLVGGFDE DGTPHLFQTD PSGIFYEWRA NTTGRSSQPV RDYMEKHADE ILTIADAAAA IKHIVRTLVS VSSLNHTQME VAVLKYRQPL RMIDHQVLAD LERTVRREIE DEAEASRRPR AP
Source	E.coli
Target Names	Prosalpha4T2
Protein Names	Recommended name: Proteasome subunit alpha type-7-1B EC= 3.4.25.1 Alternative name(s): Testis-specific alpha4-t2 proteasome subunit Testis-specific proteasome 28 kDa subunit 1B
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.