



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At4g22430 (At4g22430)

Product Code	CSB-EP638158DOA
Abbreviation	At4g22430
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1PE57
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEKPNNTITD VLEGIVTEIL VRLPLRSISR FKSVSQTWKS AIESVYFRRL FVSLHQNTSS SWSLLLRKEE FIDFHGCGTW GLPKSLGSYI QCMELDGKFE YMWFSGSNGL ILMHRKLGW KNYVGNPVLQ QWVEIPACPG SYTFFCGVVT GVDEVGVVLS FKVVKSGNMF LNKGEMYMPL YVYSSETGFW IHKEVVCPVR LPNFYDPISL NGTLYFSQRG DSYNRRPGLM VLDFYGGPKD CHFIPLPDHA LNRNKTCLAT SSGFVMYIKT LAQPGGNLLK VWRLIDDSAW QLMWEVSIPF IGCYAPMSMH PFDRNIVYLW SHDNCYLMSF NLQTQNYKIF GDESKHHDCY INHRTCEKHM YKISRPYSVS EYEGPIILQQ WVLSRWMQSV PRPPEVEMID TTSLLSLVVK NEE
Source	E.coli
Target Names	At4g22430
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At4g22430
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.