



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g62350 (At1g62350)

Product Code	CSB-MP638162DOA
Abbreviation	At1g62350
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1PFH7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSKEGLIAAK ELKRLQTQSV RLDRFIGSHV SRLKSDLVS VLAEFQRQNN VFLCMKLYEV VRREIWYRPD MFFYRDMLMM LARNKKVDET KKVWEDLKKE EVLFDQHTFG DLVRGFLDNE LPLEAMRLYG EMRESPDRPL SLPFRVILKG LVPYPELREK VKDDFLELFP GMIVYDPPED ICEDSDEEAR TDSLE
Source	Mammalian cell
Target Names	At1g62350
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At1g62350
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.