



Recombinant Meiotic recombination protein spo-11 (spo-11)

Product Code	CSB-BP638723CXY
Abbreviation	spo-11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22236
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MYEYSFNPNI DHEPGSVESQ QSTIYSDSDD SDDSFLDDEV IPPKEQAMRK IEFALADIKR QMDNKEKSLT LRISTSKSHF CLRYTAKRKG KLDRDLHCLH QVYDLLENDK RSTKRELYYE HKAUVYGNQKY LDSSIKSICE LLNESRANLN ILSCGRGIIR GAITFLVENV GVIDARVQEV LITDALLFSN IISEADFILV VEKDTTFQKL MDENFQAMFP RGILATSKGY PDIATRNVLK MLSEKRKFPI YGLFDADPHG IEIYLYKYG PTKEFAEGRG AFVPTIEWIG LFPTDFHRFT IDQSQCPLPV RTDFVKIEKM IPRSIQLGEI VVTRELDWMI QNKFKMELES INMCGQEYMA RFLIAPRVMS IEKEIPIQPE TIINEYHEDS QCSSLSTDDDR EAKDDDYIDS DAEKFKQNMI DNDS
Source	Baculovirus
Target Names	spo-11
Protein Names	Recommended name: Meiotic recombination protein spo-11
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.