



Recombinant Zinc metalloproteinase nas-21 (nas-21)

Product Code	CSB-YP638727CXY
Abbreviation	nas-21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22401
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	NKIVMR VGGSPETKRL ERSKRQALRM DNEPRWPRGT INYFFDEQRF DENSRATVLR AMEKISNHTC IKFSPKDARI KLRIVSDKGC QAAIGRVGGD QQYLSFPTSC YSVGSASELI HVIGFLHSHQ RADRDEYLKL NLQPWRLNDW FQTMQYKKYL DQWWIVPYDY GSIMQYHDS NEYGPKNSKY FRTMGSQIPS YFDYLMINEY YQCSCEGEEQ INCKNRGYPN PGNCSECNCPL LESSEEWKNI TLNLDAGYMS LQNGTKLHQI DVFVFRYLSIS APANKTIEVD IREITGVECN YGCYIGGIEV KTHEDRRMTS PRLCCKNDNE IYKSRNNPTI VMAFNSEGLD KYNIFYRFTD
Source	Yeast
Target Names	nas-21
Protein Names	Recommended name: Zinc metalloproteinase nas-21 EC= 3.4.24.- Alternative name(s): Nematode astacin 21
Expression Region	25-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.