



Recombinant Alpha-tubulin N-acetyltransferase 2 (W06B11.1)

Product Code	CSB-MP638737CXY
Abbreviation	W06B11.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q23192
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEIAFDLSTI FTDNIQRLTR TDLLKYGPKR YWAVAQSIDC LGEMSSKFHG WKRVTMYDK IVDHDEEQTT YIMWEKVNGS KSILKGLLRV GYKTLYLTDN EQNQYMEKAM CILDFFVPT EQRSGNGFKM FDEMLKAENV TVDQCAFDKP SAALQQFLEK YYDRKDLVWQ SNKYALCSNF FIGRHPTVPF TPRQTKRASR ASSAVSSHAS SRNTSPIGRN RPRHDSVADL MRQDMLAGVR AEVDPNSPTG LKNARDFGHR RIW
Source	Mammalian cell
Target Names	atat-2
Protein Names	Recommended name: Alpha-tubulin N-acetyltransferase 2 Short name= Alpha-TAT 2 Short name= TAT 2 EC= 2.3.1.108 Alternative name(s): Mec-17 like protein
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.