



Recombinant Uncharacterized WD repeat-containing protein ZC302.2 (ZC302.2), partial

Product Code	CSB-BP638739CXY
Abbreviation	ZC302.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q23256
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNPERQETIS PKNVPFQPVP TPNQQSLQSR MLESNQAVPD AFQKVAAWNH YNAVATPPIG VPQTARPPSN QSPHPNPPGY PYQSHQQQFH PLAIAFQQPF GLPQQVIYPP PPGFTPHRV QVSPVGLLGP PGFFPQFAGS TTSQGNPSEV SVRTKSAEGT TGPIAPSITT KPTSTIQVAP PRDPVAPTTT SSGITKKPEN GEFSLVKTIS GHTKSVSVIK FSYCGKYLGT GSADKQIKVW NTVDMTYLQT LASHQLGIND FSWSSNSQFI ASASDDTTVK IFDVISGACL RTMRGHTNYV FCCSFNPQSS LIASAGFDET VRVWDFKTGL CVKCI PAHSD PITSISYNHD GNTMATSSYD GCIRVWDAAS GSCLKTLVDT DHAPVTFVCF SPNGKYLLSA QLDSSLKLWD PKKAKPLKYY NGHKNKKYCL FANMSVPLGK HIISGSEDGR ILVWSIQTKQ IVQILEGHTT PVLATDSHPT LNIIASGGLE PDNVIRIWRN
Source	Baculovirus
Target Names	wdr-5.3
Protein Names	Recommended name: Uncharacterized WD repeat-containing protein ZC302.2
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.