



Recombinant Rat Granzyme B (Gzmb)

Product Code	CSB-EP010082RA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18291
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	IIGGHEAKPH SRPYMAYLQI MDEYSGSKKC GGFLIREDFV LTAAHCSGSK INVTLGAHNI KEQEKMQQII PVVKIIPHPA YNSKTISNDI MLLKLKSKAK RSSAVKPLNL PRRNVKVKPG DVCYVAGWGK LGPMGKYSdT LQEVELTVQE DQKCESYLKN YFDKANEICA GDPKIKRASR RGDSGGPLVC KKVAAGIVSY GQNDGSTPRA FTKVSTFLSW IKKTMKKS
Source	E.coli
Target Names	Gzmb
Protein Names	Recommended name: Granzyme B EC= 3.4.21.79 Alternative name(s): Fragmentin Natural killer cell protease 1 RNKP-1
Expression Region	21-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface nonself antigens, usually peptides or proteins resulting from infection by intracellular pathogens. This protein is crucial for the rapid induction of target cell apoptosis by CTL in cell-mediated immune response.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.