



Recombinant Transcription factor lin-26 (lin-26)

Product Code	CSB-EP638885CXY-B
Abbreviation	lin-26
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q27355
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLSKFVVVEV SNSNNTMTLV EYLVTHGFDK FAVVTDLVVDH PSFKYKDGSS SPESPSTTAS TAAQHTPPRT AVSTPTSINT PVPPHQNKQR HSIDKIAANL STKKVSPSSI EKQLQRTSHN PLHQLSTPHA LSQQLKLEEE QHKINQMNIQ KKEQERQQAE IQRILLQQAN AAQINHSFGL ERLTPEYDDN HHSETISKAS SEDLKTEPDS TDFGLGTSDD QVRASMLHLL HPVFAPAFGM LDAENIFGAA SKPTTPASKR RNTDSNGAPS KKHRWLPVNE LEESRSSRGK NCGRVHCKAT YKCALCGKPT TLNSTGSRWN LLRHVIMIHS DCKPYKCWDC DFTGIKSNVI SHARQCRHNA EDAHDITTDE MRAEWNLRLH ECFPDYVRAK ERGWQPEEVT VKKEEVEESP TLVKQELTLV KQEPTFAEQL EPMAQPLV
Source	E.coli
Target Names	lin-26
Protein Names	Recommended name: Transcription factor lin-26 Alternative name(s): Abnormal cell lineage protein 26
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>