



Recombinant Uncharacterized ORF1 protein (ORF1)

Product Code	CSB-BP638911NAAM
Abbreviation	ORF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27YH0
Product Type	Recombinant Protein
Immunogen Species	Nora virus
Purity	≥85% (SDS-PAGE)
Sequence	MINNQTNKKG PQLERVHFGS TQVVGKSTKR RQRGTKLDIE YTVRRNDAPK EQKFLISEIF DEKLDKQIKH EKKQNHTFIK PKLSLVIKEE QHITKKVLRG KERAATHAFM KEMVESNKIQ PSWNVEYEKE IDEVDLFFMK KKTkPFSGFS IKELRDSLIV QSDDKNMAQP TVMSSIDEIV TPREEISVSA ISEQLASLME RVDKLEKMNA ALEEENKQLK KERATIKSV KKEAKKIKQE KPQIVKKTQH KSLGVNLKIT KTKVVGQEQC LEIENTQHKK FVEKPSMPLK VSKKMTEHQL KKTIRTWYEF DPSKLVQHQK EVLNSVVTNT TFADKVRETG IPKQKIRYVA KPPAEEKRSI HFYGYKPKGI PNKVVWNWVT TGTAMDAYEK ADRYLYHQFK REMMIYRNKW VKFSKEFNPY LSKPKMVWEE NTWEYKYKTD VPYNFILKWR QLVQTYKPNT PIQADWYKIS QKQQC
Source	Baculovirus
Target Names	ORF1
Protein Names	Recommended name: Uncharacterized ORF1 protein
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.